

Computed Tomography
SOMATOM Force



Is CT nog het stralingskanon van de radiologie?



SOMATOM Force

Key components

Vectron™ X-ray tubes

- 2 x 1,300 mA @ 70 - 150 kV
- 2 x 120kW generator



Stellar^{Infinity} detectors

- 2 x 96 rows (2 x 6 cm)
- 3,120 detector channels



High-speed patient table

- 737 mm/s speed
- up to 22 cm perfusion range
- up to 80 cm dynamic range
- Turbo Flash or High Pitch scan



Low-dose early detection with Tin Filter CT imaging at the dose of conventional X-ray¹⁾

SOMATOM Force

Collimation: 192 x 0.6 mm

Scan time: 0.4 s

Scan length: 200 mm

RoI

Sn1

CTI

DLI

Eff.

Lur

imp

allow doses of conventional X-ray¹⁾



Table 1

Radiology 2008;248:254-63

Adult Effective Doses for Various Diagnostic Radiology Procedures

Examination	Average Effective Dose (mSv)	Values Reported in Literature (mSv)
Skull	0.1	0.03–0.22
Cervical spine	0.2	0.07–0.3
Thoracic spine	1.0	0.6–1.4
Lumbar spine	1.5	0.5–1.8
Posteroanterior and lateral study of chest	0.1	0.05–0.24
Posteroanterior study of chest	0.02	0.007–0.050
Mammography	0.4	0.10–0.60
Abdomen	0.7	0.04–1.1
Pelvis	0.6	0.2–1.2
Hip	0.7	0.18–2.71
Shoulder	0.01	...
Knee	0.005	...
Other extremities	0.001	0.0002–0.1
Dual x-ray absorptiometry (without CT)	0.001	0.001–0.035
Dual x-ray absorptiometry (with CT)	0.04	0.003–0.06
Intravenous urography	3	0.7–3.7
Upper gastrointestinal series	6*	1.5–12
Small-bowel series	5	3.0–7.8
Barium enema	8*	2.0–18.0
Endoscopic retrograde cholangiopancreatography	4.0	...




Early detection with low-dose lung scans

Outstanding image quality even in obese patients

SOMATOM Force

Collimation: 192 x 0.6 mm
Scan time: 2.0 s
Scan length: 274 mm
Rotation time: 0.5 s
110 kV, 79 mAs
CTDI_{vol}: 4.13 mGy
DLP: 134 mGy cm
Eff. dose: 0.28 mSv

Early detection low-dose lung scan in an **obese patient** with a **very low dose of 0.28 mSv** with outstanding image quality.



Cardiac imaging with high heart rates – ultra-low-dose Turbo Flash Spiral at a heart rate of 117 bpm

SOMATOM Force

Collimation: 2 x 192 x 0.6 mm

Scan time: 0.12 s

Scan length: 93 mm

Rotation time: 0.25 s

80 kV, 289 mAs

CTDI_{vol}: 1.33 mGy

DLP: 19.7 mGy cm


Eff. dose: 0.28 mSv

HR: 117 bpm

CM: 30 mL

Rule out of coronary disease

within 0.12 s at average heart rate
of **117 bpm** with **Turbo Flash mode**.




Is CT nog het stralingskanon
van de radiologie?

**Dosis reducerende
technologie in de
Siemens SOMATOM Force**




CARE technology

Combine Applications Reduce Exposure



More precise configuration of dose modulation

- Configuration for different body habitus and organs individually
- New and more specific modulation curves
- Definition of threshold values for Dose Alert according to IEC regulations




CARE Dose4D and CARE kV



Dose relevant scan parameters


- Patient's body habitus
- Effective current (mAs)
- Voltage dep. on examination (kV)

Current modulation w/ CARE Dose4D



- Automatic patient and organ specific protocol adaption
- Automatic modulation in X,Y & Z direction in real time



Voltage setting with CARE kV



Example: Abdomen @ 100 kV




Dose (CTDIvol):	9.7 mGy
Dose reduction compared to:	
• Siemens 120 kV stand. prot.:	14 mGy
• Competitor's 120 kV protocol*:	25 mGy




*Source: Internal data evaluation based on anonymous assessment on SRS connected scanners. Clinical image courtesy C. McCullough, Mayo Clinic, Rochester, MN, USA * Sagara Y et al. Abdominal CT: comparison of low-dose CT with adaptive statistical iterative reconstruction and routine-dose CT with filtered back projection in 53 patients. AJR Am J Roentgenol. 2010 Sep;195(3):713-9.

The newest generation in iterative reconstruction

Non ADMIRE



ADMIRE




ADMIRE in three steps

- **Raw data statistical modeling:**
Statistical weighting of all projections in raw data space
- **Model based noise cancelation:**
Noise cancelation based on an intelligent approach in image space
- **Advanced system modeling:**
Advanced modeling is the base for the forward projection and eliminates artifacts

ADMIRE

Excellent image quality



Patient with dissection

Low kV imaging with reduced contrast media (20 mL)

SOMATOM Force

Collimation: 2 x 192 x 0.6 mm

Scan time: 1.1 s

Scan length: 740 mm

Rotation time: 0.25 s

80 kV, 140 mAs

CTDI_{vol}: 2.09 mGy

DLP: 155 mGy cm

Eff. dose: 2.32 mSv

Gated, examination with
20 mL¹⁾ of contrast media
without breath-hold
in patient with
limited function in right kidney.




Two steps ahead in Freezing Motion

“Free-breathing” CT imaging

Ultra-high-speed imaging with the unique SOMATOM Force in combination with ADMIRE

Next generation Dual Source CT (DSCT) with

- Two Vectron™ tubes
- Two Stellar^{Infinity} detectors with 25% more detector channels (1,840)
- Up to 22 lp/cm (0.24 mm) spatial resolution
- Advanced Modeled Iterative Reconstruction (ADMIRE)¹⁾ for lowest possible dose




Up to 66 ms native temporal resolution utilizing DSCT at 250 ms rotation speed

1) Image quality as defined by low contrast detectability using a model observer method for evaluation. Equivalent low contrast detectability can be achieved with 80% to 85% less dose using ADMIRE at highest strength level for thin (0.6 mm) reconstruction slices in measured and simulated body and head phantoms for low contrast objects with different contrasts. See ADMIRE date sheet for further information. In clinical practice, the use of ADMIRE may reduce CT patient dose depending on the clinical task, patient size, anatomical location, and clinical practice. A consultation with a radiologist and a physicist should be made to determine the appropriate dose to obtain diagnostic image quality for the particular clinical task.

“Free-breathing” CT imaging

No anesthesia, no sedation¹⁾ for babies/kids



Baseline (previous high-end CT)
Standard low-dose thorax CT
3 mSv effective dose
11 s scan time



Follow-up (SOMATOM Force)
Turbo Flash ultra-low-dose thorax CT
0.4 mSv effective dose; **0.3 s** scan time
Pediatric imaging without breath-hold and sedation!¹⁾

Courtesy of UMM, Mannheim, Germany

¹⁾ The inherent temporal resolution – the ‘native’ temporal resolution acquired by the scanner – is highly important to freeze patient motion, e.g. in lung exams or in patients who cannot hold their breath long enough. This is also important, in pediatric CT where it also can help reducing the need for potentially harmful sedation.

Preventive care in cardiac imaging

Ultra-low-dose Turbo Flash Spiral with 0.09 mSv

SOMATOM Force

Collimation: 2 x 192 x 0.6 mm

Scan time: 0.2 s

Scan length: 142 mm

Rotation time: 0.25 s

70 kV, 128 mAs

CTDI_{vol}: 0.36 mGy


DLP: 6.6 mGy cm

Eff. dose: 0.09 mSv

HR: 56 bpm

BMI: 24

Rule out of coronary disease **within 0.2 s** and with ultra-low-dose with **only 0.09 mSv.**



Two steps ahead in Decision Making


Precise Dual Energy quantification

Dual Energy (DE) imaging with improved energy separation

Vectron™ tube plus Stellar^{Infinity} detectors

- New energy pairings, e.g. 90 and 150 kV Sn for imaging of obese patients
- Up to 35 cm DE field of view (FoV)
- Up to 258 mm/s DE scan speed
- Tin Filter (Selective Photon Shield II) utilization

Up to 30% increased energy separation for better DE imaging outcomes



Precise Dual Energy quantification – pancreatic head tumor illustrated with different DE reconstructions

SOMATOM Force

Collimation: 2 x 128 x 0.6 mm

Scan time: 4.35 s

Scan length: 207mm

Rotation time: 0.5 s


100/Sn150 kV, 117/60 mAs

CTDI_{vol}: 6.88 mGy


DLP: 162.9 mGy cm

Eff. dose: 2.4 mSv


Advanced diagnostic information:
Monoenergetic Plus allows energy spectrum shift to lower levels which **significantly improves** image quality and **pancreas-to-lesion contrast¹⁾**.



Mono+ 55 keV



Mixed



VNC

Courtesy of MUSC Medical Center, Charleston, USA


1) Frellesen et al. (2015): Dual Energy CT of the pancreas: improved carcinoma-to-pancreas contrast with a noise-optimized monoenergetic reconstruction algorithm.

Precise DE quantification – tissue differentiation and therapy response assessment of liver cancer

SOMATOM Force

Collimation: 2 x 192 x 0.6 mm
Scan time: 4 s
Scan length: 182 mm
Rotation time: 0.5 s
100/Sn150 kV, 202/190 mAs
CTDI_{vol}: 11.99 mGy
DLP: 254 mGy cm
Eff. dose: 3.81 mSv

Dose neutral Dual Energy enables advanced diagnosis of liver metastases and cyst tissue.



Kidney-friendly scanning – change treatment decisions with more comprehensive 4D imaging

SOMATOM Force

Collimation: 48 x 1.2 mm

Scan time: 36 s

Scan length: 222 mm

Rotation time: 0.25 s

70 kV, 200 mAs


CTDI_{vol}: 43.46 mGy

DLP: 905 mGy cm

Eff. dose: 13.6 mSv

CM: 12 mL

Kidney-friendly dynamic angiography of EVAR, Type II endoleak scanned **with only 12 mL¹⁾** of contrast media in 72 year old patient.



SOMATOM Force

Two steps ahead

Freezing motion – avoiding artifacts
Dual Source CT with Turbo Flash Spiral

Kidney-friendly scanning
Vectron™ X-ray tubes with maximum power at low kV

Perfusion imaging in routine
Dynamic imaging with low kV and unique Adaptive Dose Shield



Outstanding quantitative certainty
Fast Dual Energy acquisition with improved spectral separation

More channels. More coverage. Excellent image quality.
Stellar^{Infinity} detectors – with anti-scatter 3D collimator grid

Low-dose early detection
Tin Filters with ADMIRE

Slices	2 x 192 (2 x 576 recon)	mA	1,300 mA@70 kV 1,300 mA@80 kV	Rotational coverage	184 mm/rot	kV	70 - 150 kV 10 kV Steps
Max. scan speed	737 mm/s	Temp. res.	66 ms	Power	240 kW (2 x 120 kW)	Channels	1,840

Is CT nog het stralingskanon van de radiologie?



Dank voor uw aandacht

Computed Tomography
Listen.
Care.
Improve.